

## *EN General notes on alcohol testing - Instructions for use*

The alcohol test is designed for single use and is used to determine the alcohol content in the exhaled air. This test allows you to detect blood alcohol with an accuracy of 0.2‰, or 0.1 mg/L in exhaled air. The alcohol level tends to increase up to about an hour after consumption. Factors such as individual physical condition, rate of consumption, type of alcohol, and food intake can affect alcohol level.

### **Attention!**

The results of this alcohol test cannot be used as evidence in legal proceedings, in particular to determine the state of intoxication. You are solely responsible for interpreting the test results. These results cannot be used to challenge the results of tests performed by or at their request. The manufacturer assumes no responsibility for use of the product that does not comply with the instructions.

The manufacturer and distributor do not accept any legal responsibility for the outcome of the test.

To obtain the most accurate result possible, the measurement should be taken at least 15 minutes after the last alcohol consumption. Smoking and consuming certain foods and beverages can alter the outcome; even in these cases it is advisable to wait at least 15 minutes. In any case, the rule must be:

**IF YOU DRINK YOU DON'T HAVE TO DRIVE ! IF YOU DRIVE YOU DON'T HAVE TO DRINK!**

### *Components*

The alcohol test consists of a glass tube, adsorption granules, filters, aluminum foil and two protective caps.

### **Attention!**

The components of the alcohol test can be dangerous if used improperly. The glass tube may break and cause injury, other parts may be swallowed. After use, the test should be disposed of in household waste. The test used does not pose an environmental hazard according to EU directives. It is recommended that you do not leave used tests in places that are accessible to unauthorized persons.

### *Instructions:*

1. Press the plastic caps on both sides of the alcohol test until you feel resistance to pierce the film

Protective.

2. Inhale deeply and blow air into the mouthpiece following the arrows on the tester (twice for about 10 seconds each). Make sure that the air comes out the other side of the tester.

3. Read the test result after 2 minutes, but no later than 5 minutes. 4. If the color in the middle of the glass tube changes, it means that there is alcohol in the blood.

5. Compare the test color with the scale on the packaging. The most precise results are obtained in daylight. A similar value of 0.2‰ indicates a blood alcohol concentration of 0.2‰ or more.

An alcohol level of 0.2‰ can impair driving ability and poses a danger to one's life and health, as well as a risk to other road users.

**Attention!**

Please comply with the legal regulations regarding the maximum blood alcohol limits in your country!

*Precautions:*

1. Keep the alcohol test out of the reach of children.
2. It is forbidden to remove the measuring substance from the alcohol test.
3. Do not take the measurement immediately after using a glass cleaner containing alcohol, as the vapors can alter the result.
4. Avoid direct contact of the measuring medium with eyes and skin.
5. The measuring substance must not be inhaled or swallowed.
6. Do not dispose of the alcohol test in a fire.
7. If the measuring substance comes into contact with the eyes, open the eyelid well and rinse thoroughly with water. Seek medical attention immediately, an eye examination is advisable.
8. In case of skin contact, rinse the affected area thoroughly with water. In case of skin irritation (redness, rashes, allergic reactions, wounds), consult a doctor and dermatologist immediately.
9. If swallowing the measuring substance, drink plenty of water, but do not induce vomiting. Seek medical attention immediately.
10. If you inhale the measuring substance, go to fresh air immediately. If necessary, consult a doctor.

Storage tips: 1. Protect from sunlight and high temperatures.

2. Store in a dry and safe place out of the reach of children.

**MADE IN EU**

***Authorized Distributor:***

My Future SAS Innovation

37- 41 Boulevard Dubouchage

06000 Nice, France

E-mail: [info@myfutureinnovation.com](mailto:info@myfutureinnovation.com)

The alcohol test is produced in Europe.